
ERKOFLEX-95

DIN EN ISO
13485: 2003

Erkoform-3/3d
Erkopress 300 Tp/-ci

1,5 mm, 155 °C
2,5 mm, 140 °C
4,0 mm, 130 °C

1,5 mm, 1:50 min
2,5 mm, 2:20 min
4,0 mm, 3:30 min

Erkoform-RVE
Erkopress 200 E

1,5 mm, 1:15 min
2,5 mm, 1:55 min
4,0 mm, 2:45 min

Lagerung • storage • stockage

ERKODENT Erich Kopp GmbH

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72285 Pfalzgrafenweiler, Germany

Tel. + 49 (0) 74 45/85 01-0

Fax + 49 (0) 74 45/20 92

info@erkodent.com • www.erkodent.com

DE ERKOFLEX-95

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EN ERKOFLEX-95

Chemical characteristics:

EVA: ethyl-vinyl-acetate
PE, insulating foil: polyethylene
Harmless to health, BgVV (Federal Health Office)
listed thermoplast with tested biocompatibility.
Insoluble in water, inactive, harmless to ground
water.

Disposal/Recycling:

Household waste (combustible, low in pollutants and
high in energy)/general recycling for plastic.

Field of use:

(see brochure Thermoforming technique and
material card) Resilient, durable, highly comfortable
occlusal splints which can be bonded with the
Erkoflex-95 material in sticks (use the fusing gun) or
occlusal splints adjustable by occlusal warming.

Material properties:

Flexible material with high elasticity and dimensional
stability. Can be adjusted and bonded with the fusing
gun (Erkoflex-95-sticks). Does not bond to
acrylate. Hardness app. 95 Shore A.

Finishing: (see brochure Thermoforming technique)
Finishing set Quick 2 (110 877) or scissors. The
spacer foil will be taken off only after thermoforming
and finishing.

Plastification: (for machines without heating time
control) Place the foil in the devices in such a way
that the spacer/insulating foil is pointing towards the
model. Heat until foil sags slightly and becomes
clear, then thermoform.

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